

VILLA LABECO s.r.o. Chrapčiakova 1 052 01 Spišská Nová Ves

ISOTACHOPHORESIS

APPLICATION NOTE No. 11

DETERMINATION of SORBIC ACID in FOOD

MAIN FEATURES:

Sorbic acid as a preservative can be determined quantitatively by isotachophoretic separation in drinks, wines, jams, etc. Little or no sample pre-treatment, and high reproducibility of the results are typical features of this determination. Simultaneously it is possible to analyse other anions: benzoic, formic, lactic, PO₄, citric, etc.



Fig. 1 : Isotachophoreogram of syrup concentrate for preparation of lemonade BELLA (fy White Lady) A – Isotachophoreogram from conductivity detector of preseparation column. Citric acid was quantitatively determined. Only the zones in the box were analysed in analytical column.

- B Isotachophoreogram from conductivity detector of analytical column
- C Isotachophoreogram from UV-detector. 500 mg/kg is determined concentration of sorbic acid
- D Isotachophoreogram from UV-detector, syrup was free of sorbic acid.

terminating electrolyte (TE): 5.10⁻³ M caproid acid

 I_1 = 250 $\mu A, I_2$ = 50 $\mu A,$ V=30 $\mu I,$ sample dilution 1 :200



- Fig. 2 : Isotachophoreogram of wine with adding 100 mg/l of sorbic acid A – record from conductivity detector B – record from UV-detector (254 nm)
- Conditions : as in Fig.1., only leading electrolyte in analytical column was with half concentration for increasing sensibility.

Leading electrolyte : 5.10^3 M Cl + 2,5.10⁻² M β -alanine + 0,1% HEC+ pH=3,9 I_2 = 30 μ A; sample dilution 1 : 25

CZE and ITP analysers are produced by : Villa Labeco s.r.o., Chrapčiakova 1, 052 01 Spišská Nová Ves, Slovakia www.villalabeco.sk